

# Non-Invasive Phase-Based Neurostimulation for Enhancing Attention

## Background: Rhythms of Attention and Phase Entrainment

Human attention is closely tied to brain oscillations – rhythmic patterns of neural activity at specific frequencies <sup>1</sup> <sup>2</sup>. For example, theta-band (~4–7 Hz) activity in fronto-parietal networks is linked to executive attention and working memory, alpha-band (~8–12 Hz) oscillations in parietal and frontal regions correlate with selective attention and inhibition of distractions, and gamma-band (>30 Hz) oscillations are associated with active information processing <sup>2</sup>. Crucially, *phase* relationships (the timing alignment of these oscillations across brain regions) determine how effectively different areas communicate. Research shows that when distant regions oscillate in-phase (at ~0° phase lag), information transfer and cognitive performance can improve, whereas out-of-phase activity (e.g. 180° lag) can hinder it <sup>3</sup>. These insights have spurred development of **phase-based neurostimulation** techniques – using external stimuli to **entrain** brain rhythms at particular frequencies and phases – as potential tools to boost attention. We focus on non-invasive methods, including **transcranial alternating current stimulation (tACS)**, **rhythmic sensory entrainment** (through visual or auditory rhythmic stimuli), and **EEG neurofeedback**, highlighting human studies and protocol parameters (frequency, phase offset, timing, etc.) relevant to an “attention phase trainer” tool. (Invasive findings, such as intracranial recordings of phase-specific effects, are mentioned only to inform mechanisms, not as direct interventions.)

## Transcranial Alternating Current Stimulation (tACS) for Attention

tACS involves applying weak sinusoidal currents to the scalp at specific frequencies to rhythmically **entrain endogenous brain oscillations** <sup>4</sup>. By aligning neural firing to an external oscillation, tACS can modulate the phase and power of brain rhythms linked to attention. Multiple studies show frequency- and phase-specific effects:

- **Theta-frequency tACS (~6 Hz)** over frontoparietal networks can enhance cognitive performance when the stimulation of two regions is phase-synchronized. In a landmark study, Polanía et al. applied 6 Hz tACS to left frontal and parietal cortices either in-phase (0° phase lag) or out-of-phase (180°) <sup>5</sup>. Only the **in-phase (synchronized) theta tACS** significantly **improved reaction times** in a visual working memory task, whereas anti-phase stimulation **worsened performance** <sup>3</sup>. This demonstrated a causal role of theta phase-coupling between frontal and parietal areas in enhancing cognitive function. An attention-training device could leverage this by driving theta oscillations coherently across frontal and parietal sites (e.g. F3/P3 or F4/P4 montages) to strengthen frontoparietal network connectivity. Typical tACS protocols use **20-minute sessions** at comfortable intensities (~1–2 mA), often repeated over multiple sessions for lasting effects <sup>6</sup>. Notably, **phase offsets** can be tuned: a 0° relative phase is hypothesized to facilitate communication between regions, whereas 180° might simulate distractor conditions or could be avoided in pro-attentional training <sup>7</sup>.

- **Alpha-frequency tACS (~10–12 Hz)** has been used to target the brain's attentional gating mechanisms. Alpha oscillations are thought to reflect top-down suppression of irrelevant inputs <sup>8</sup> . For instance, **10 Hz tACS over frontal eye fields (FEF)** was shown to **mimic an endogenous attention state**, altering visual processing <sup>9</sup> <sup>10</sup> . In a recent study, intermittent alpha-tACS over bilateral FEF changed participants' visual contrast perception: it **improved sensitivity to faint stimuli** while slightly **impairing perception of very high-contrast stimuli**, effectively altering the psychometric curve for visual detection <sup>11</sup> <sup>12</sup> . The authors propose that alpha tACS **"strengthened a fronto-occipito-parietal network in the alpha band,"** mimicking top-down attention and enhancing phasic gating in visual cortex <sup>10</sup> . This suggests that an attention trainer could use alpha-tACS over frontoparietal sites to help users better filter out distractions (by reinforcing neural inhibition patterns similar to focused attention). Indeed, another study found that 10 Hz occipital tACS entrained ongoing alpha rhythms and could induce periodic fluctuations in visual attention performance <sup>13</sup> (e.g. aligning with alpha phase could predict moments of impaired target detection <sup>13</sup> ). Care must be taken, however, as alpha-tACS might also reduce temporal resolution of perception (slowing the rate of sampling) <sup>14</sup> , so its use may be optimal when the goal is to stabilize attention and block out rapid distractions.

- **Gamma-frequency tACS (~30–50 Hz)** is associated with enhancing local processing and arousal. Gamma tACS has shown promise in improving certain attention-related metrics, especially in clinical populations. A notable example applied **40 Hz tACS over the right frontoparietal (FP) network** (electrodes at F4 and P4) in a series of 12 sessions (20 minutes each at 2 mA) for individuals with attention deficits <sup>6</sup> . This protocol led to **significant improvements on a sustained attention test (TOVA)** – specifically, variability of reaction time decreased, signal detection (d') improved, and commission errors (false alarms) dropped <sup>15</sup> . Mean reaction time and omission errors remained stable, indicating that the intervention mainly improved consistency and impulsivity control rather than speed <sup>15</sup> . The authors conclude that **40 Hz tACS over the right FP network can improve attention, particularly in terms of more consistent performance and better inhibitory control** <sup>15</sup> . This aligns with the idea that gamma oscillations may boost the "signal" of task-relevant activity and executive control. For a product, gamma tACS could be delivered either globally (e.g. frontoparietal as above) or more focally (e.g. right frontal cortex) to promote alertness and executive attention. It's worth noting that gamma tACS effects can be subtle; one study in healthy adults found no significant working memory benefit from 40 Hz stimulation alone <sup>16</sup> , suggesting that baseline cognitive status or concurrent tasks might modulate its efficacy. Nonetheless, gamma entrainment remains a key feature to consider, especially given reports that **40 Hz rhythmic stimulation can engage widespread networks** (indeed 40 Hz visual flicker has even been explored for cognitive enhancement in Alzheimer's disease models, though that is beyond our scope).

**Personalization and Timing:** A recurring theme is the **individual variability** in optimal stimulation parameters. Recent work suggests that each person may have a "preferred" frequency for entrainment – essentially a resonance frequency of their attention networks <sup>17</sup> . In a study combining tACS with auditory tasks, when no strong external rhythm was present, participants showed maximal behavioral entrainment at different frequencies centered around the theta/alpha range, varying individually <sup>18</sup> . This implies that an attention-training device might benefit from a calibration phase: e.g. measuring a user's EEG to find their dominant alpha or theta frequency, and then tuning the tACS to that frequency for maximal effect <sup>17</sup> . Timing the stimulation relative to cognitive events is also critical. For instance, the FEF alpha-tACS study delivered pulses in an **intermittent** manner, avoiding stimulation during the exact moment of target stimuli to prevent masking or sensory side effects <sup>19</sup> <sup>20</sup> . Future "phase trainer" tools could incorporate **closed-**

**loop timing**, where stimulation is triggered at specific phases of the user's ongoing brain oscillation or at specific task intervals, to optimally synchronize with natural attentional cycles. Such closed-loop tACS approaches are still experimental but conceptually could amplify the benefits by ensuring the external rhythm aligns with the brain's current phase state.

It is important to mention that **concurrent sensory input can compete with tACS-driven entrainment**. A recent **2025 study (Cabral-Calderin & Henry)** found that when people received **rhythmic auditory stimuli together with tACS**, the **sensory rhythm dominated the behavioral outcomes**, effectively overshadowing tACS unless the sensory input was very subtle <sup>21</sup> <sup>22</sup>. Only when the auditory rhythm was removed or had no informative pattern did tACS alone significantly entrain behavior, and even then the best tACS frequency differed across individuals <sup>23</sup> <sup>17</sup>. This finding highlights that an attention training device should consider **modality interactions**: if using tACS, one might train in a quiet or non-flickering environment (to let the electrical entrainment take hold), or conversely if combining with rhythmic sensory cues (see next section), design the protocol such that tACS and sensory rhythms complement rather than conflict. In sum, tACS offers a versatile way to externally drive brain oscillations relevant to attention – with the ability to target specific frequencies (theta/alpha/gamma) and phase relationships (synchronized vs. desynchronized) between regions. Effective use will likely involve **personalized frequency tuning, appropriate montage selection (e.g. fronto-parietal electrode placements to engage dorsal attention networks), and careful timing to integrate with cognitive activities**.

## Rhythmic Sensory Entrainment (Visual and Auditory)

Another non-invasive route to influence brain oscillations is through **rhythmic sensory stimulation**, such as flickering visual stimuli or patterned sounds that naturally induce brainwave entrainment. When the brain is exposed to a repetitive stimulus, neurons tend to synchronize their firing to that rhythm – producing steady-state evoked potentials (SSEPs) at the stimulation frequency <sup>24</sup> <sup>25</sup>. Unlike electrical stimulation, which directly injects currents, sensory entrainment works via normal sensory pathways but can still drive oscillatory brain activity and in turn affect attention and cognition.

- **Visual Entrainment (SSVEP-based)**: Presenting flashing images or lights at specific frequencies can lock cortical activity to that frequency. This is the basis of SSVEP-based brain-computer interfaces and has been widely considered neutral in terms of cognition. However, recent evidence shows **frequency-specific visual flicker can causally alter attention performance** <sup>26</sup> <sup>13</sup>. For instance, exposing participants to **different flicker rates (e.g. 10 Hz vs 40 Hz)** in the background of a visual attention task changed how they allocated attention between left and right visual fields <sup>27</sup>. In one study, **40 Hz visual flicker led to a bias favoring rightward attention and enhanced neural markers of attentional suppression (the Pd component), compared to 10 Hz flicker** <sup>28</sup>. Meanwhile, **10 Hz flicker (alpha-band) has been shown to entrain occipital alpha oscillations and cause periodic fluctuations in visual detection performance** <sup>13</sup>. Specifically, an alpha rhythm in the visual input can impose moments of high excitability and low excitability, thereby creating an endogenous-like “sampling” cycle <sup>13</sup>. Gulbinaite et al. (2017) reported that a 10 Hz background flicker can actually **impair target detection** in a visuospatial attention task, presumably by synchronizing with and amplifying the brain's intrinsic alpha inhibition cycles <sup>29</sup>. On the other hand, properly tuned higher-frequency flicker might avoid that inhibitory regime and instead enhance processing (as suggested by the 40 Hz benefits above). These results imply that an attention training product could use **rhythmic visual stimuli to train the brain's attentional rhythm**. For example, a **40 Hz LED flicker** (in goggles or a screen) during training exercises might

stimulate gamma-band activity in visual and parietal cortex, potentially improving aspects of attention or spatial bias. Conversely, a **theta or alpha rhythmic visual cue** could be used in protocols aiming to train anticipatory attention in time – e.g. flashing at 5 Hz to encourage a user’s attentional focus to oscillate in sync, thereby learning to optimize detection at specific phases (some researchers have explored presenting cues rhythmically so that the brain “tunes” to expect events at that pace).

Any use of visual entrainment must account for **safety and comfort**. Certain frequencies in the beta range (15–30 Hz) and notably around 10 Hz are known to trigger photic epilepsy in susceptible individuals; thus a device would need controls to avoid unsafe flicker (e.g. perhaps using lower-contrast stimuli or excluding users with photosensitive epilepsy). Interestingly, **intermittent alpha stimulation** (turning flicker on and off periodically rather than continuous) has been found to produce after-effects on EEG – one study showed that intermittent 10 Hz visual flicker increased resting alpha power even after the stimulation ceased <sup>30</sup>. This suggests a strategy where a trainer tool might use bursts of flicker to not only entrain in the moment but also gradually strengthen the user’s endogenous alpha rhythm (if the goal is, say, to help an overactive brain learn to idle at a higher alpha power associated with calm focus).

- **Auditory Entrainment:** Rhythmic sound can also synchronize neural activity, particularly in the auditory cortex and beyond. **Auditory steady-state stimulation** can be done via amplitude-modulated noise or click trains at a target frequency. For attention, auditory rhythms might be used to influence alertness and temporal attention (for example, a series of tones at theta rate could encourage a listener’s attentional focus to wax and wane in that rhythm). A recent study tested combined auditory and tACS entrainment: participants listened to amplitude-modulated noise at various depths while receiving concurrent tACS, and had to detect gaps in the noise <sup>31</sup>. It was found that **behavioral entrainment to the auditory rhythm was strong and depended on the sound’s modulation depth, whereas entrainment to tACS was only detectable when the auditory rhythm was absent or very weak** <sup>21</sup> <sup>22</sup>. In the absence of a rhythmic sound, about 70% of participants could be behaviorally entrained by tACS alone, but interestingly **the optimal entrainment frequency differed between individuals** <sup>23</sup> <sup>18</sup>. When the rhythmic sound was present, it **“prevailed” in determining behavior, underscoring the potency of sensory entrainment** <sup>32</sup>. For an attention-training tool, this indicates that **auditory cues can be a powerful way to drive oscillations** related to attention. For example, a **binaural beat or isochronic tone at 4–7 Hz** might engage frontal theta networks (some studies have reported enhanced memory encoding with theta-range auditory beats). Likewise, a **40 Hz auditory click train** has been explored for stimulating gamma oscillations; in mouse models of Alzheimer’s, 40 Hz auditory stimulation induced gamma entrainment and cognitive benefits, and human trials are examining similar stimuli for cognitive support. Although evidence in healthy attentional tasks is still emerging, a plausible feature is incorporating rhythmic auditory feedback when visual flicker is impractical (e.g. if the user’s eyes need to be closed or during meditation-like attention training, gentle auditory beats might entrain brainwaves). Auditory entrainment is generally more comfortable over long periods than bright flicker, but it may be limited in spatial specificity (sound entrains more globally whereas visual flicker can be localized in the visual field).
- **Tactile Entrainment:** Though less studied, one can also use rhythmic tactile stimulation (e.g. vibration on the wrist at a certain frequency) to drive brain rhythms. This primarily engages somatosensory cortex oscillations. If relevant, an attention trainer could include a vibrating element (for example, a rhythmic tapping on a finger at alpha frequency to induce somatosensory alpha). One application might be for tasks requiring sustained somatosensory attention or to combine with

motor attention training. However, evidence for tactile-driven improvements in attention is sparse, so this would be more experimental.

In designing the sensory entrainment component of a product, **frequency selection** is key. The literature suggests **10 Hz (alpha) and 40 Hz (gamma)** are “classic” frequencies with known neural impact (alpha-flicker for entrainment/inhibition, gamma-flicker for activation). Indeed, a 2023 study explicitly compared **10 Hz vs 15 Hz vs 40 Hz visual flicker** and found each had distinct effects on attention asymmetry and EEG markers <sup>27</sup>. This frequency specificity means a device might offer several modes (e.g. “*Alpha entrainment mode*” vs “*Gamma entrainment mode*”), depending on whether the goal is to train focus and inhibition (alpha) or to stimulate processing speed and alertness (gamma). **Phase considerations** are relevant even in sensory entrainment: if combining stimuli (say, a flickering light and a beeping sound), their phase alignment could matter (synchronous audio-visual beats vs alternating). Moreover, if sensory entrainment is combined with tACS, one could try to **phase-align the tACS with the sensory rhythm** to see if that yields additive effects (for instance, delivering tACS in phase with a visual flicker might reinforce the neural response). To date, the general finding is that **sensory rhythms tend to dominate** <sup>32</sup>, but careful phasing might still prove beneficial, and future research could explore closed-loop adjustment where the phase of a flicker is continuously adjusted to maintain a desired brain-wave phase (detected via EEG).

In summary, rhythmic sensory entrainment offers a **non-electrical, user-friendly way to engage brain oscillations**. An “attention phase trainer” could use **LED goggles, screen flicker, or rhythmic auditory cues** as part of training sessions to entrain the brain at specific frequencies. The parameters gleaned from research include using moderate flicker rates (theta, alpha, gamma bands) for 5–20 minutes at a time, possibly interleaving different frequencies per session. The **benefits** can be task-specific – e.g. 40 Hz visual flicker might enhance certain visuospatial attention biases <sup>28</sup>, whereas 10 Hz might train the user to experience and overcome periodic lapses in attention. All sensory entrainment should ideally be paired with an actual attention task (e.g. the user is trying to concentrate on a game or stimuli while the flicker occurs) so that the brain learns to maintain performance under those rhythmic conditions, thereby *entraining not just the brainwaves but the cognitive skill* of syncing with or resisting the rhythm. This pairing of cognitive training with rhythmic stimulation is supported by the principle that entrainment alone can modulate brain state, but **entrainment + practice** can lead to functional improvements that outlast the stimulation period.

## EEG Neurofeedback for Attention Regulation

Neurofeedback (NFB) is a training method where users learn to modulate their own brain activity through real-time feedback, typically via EEG. Unlike tACS or sensory entrainment which drive oscillations externally, neurofeedback empowers the individual to internally adjust their brain rhythms. This approach has been extensively applied to attention-related conditions, especially ADHD, and can be considered for an attention trainer tool either as a standalone mode or combined with stimulation.

**How it works:** In a neurofeedback session, EEG electrodes monitor certain brainwave features (power in a frequency band, or ratios, or coherence between regions), and the user is presented with a feedback signal (visual animations, auditory tones, game mechanics) that reflects their real-time brain activity. Through practice, users can learn to intentionally alter that activity to reach certain targets – effectively **training their brain to sustain a more “optimal” oscillatory state for attention**.

- **Targeted Frequency Bands in Attention:** Traditional EEG-NFB protocols for improving attention often focus on the **theta and beta bands**. Individuals with attentional problems (like ADHD)

frequently show **excess theta (4–8 Hz) and/or deficient beta (13–20 Hz)** during resting states, leading to a high theta/beta ratio (TBR) associated with inattention <sup>33</sup>. A common protocol is **theta/beta ratio training**, where the goal is to **reduce theta power while increasing beta power**. For example, a child might see a rocket ship that only flies when their beta activity goes up and theta goes down; over sessions, they learn (implicitly) to produce more beta (alert, focused state) and suppress theta (drowsy, mind-wandering state). Meta-analyses have reported that such **EEG neurofeedback can improve ADHD symptoms of inattention and impulsivity** to a moderate degree <sup>34</sup>, though some studies debate the magnitude of specific vs. placebo effects <sup>35</sup>. One systematic review in 2021 found that **EEG-NFB led to significant improvements in attention and hyperactivity ratings** compared to non-active controls <sup>34</sup>, while other controlled trials suggest that when properly blinded, effects might be smaller. Nonetheless, **the ability of neurofeedback to induce lasting neural changes is well documented**, as users essentially practice the self-regulation of brain oscillations.

- **SMR and Alpha Training:** Another protocol used in attentional training is **Sensorimotor Rhythm (SMR) up-training**. SMR is a mid-range alpha frequency (~12–15 Hz) observed over the sensorimotor cortex; increasing SMR has been linked to relaxed but focused mental states. Studies have shown that training individuals to raise SMR can improve sustained attention and even aspects of motor inhibition. Similarly, **alpha neurofeedback** (usually in occipital or frontal sites) can be used either to increase alpha (for calming, as in anxiety reduction which indirectly may aid attention) or to decrease alpha (to reduce excessive cortical idling in someone who is underaroused). The choice depends on the individual's EEG profile: e.g., adults with ADHD sometimes have *excess slow theta* and *maybe excess high-beta (hyperarousal)*, so protocols can be adjusted (some protocols train down high beta to reduce restlessness, combined with training up SMR or beta in low-arousal cases). Because an attention-phase trainer device would target a broad user base, an **adaptive neurofeedback system** could first perform an EEG assessment to identify whether a user's baseline attention-related rhythms are deviant (e.g. high theta) and then personalize the feedback protocol (target that metric).
- **Phase/Coherence Neurofeedback:** Traditional neurofeedback focuses on power amplitudes, but it is also possible to feed back measures of **connectivity** (like coherence or phase locking between two electrodes). Though more technically complex, this directly ties into the “phase” aspect of attention networks. For instance, one could attempt to train a user to increase the coherence (phase synchronization) between frontal and parietal EEG channels in the theta band – effectively teaching them to internally synchronize those regions, akin to what tACS in-phase did externally. There have been exploratory studies where coherence NF was used in autism or mood disorders, but for attention this is still an emerging idea. However, conceptually, if research like Polanía's shows in-phase frontoparietal theta is beneficial, a neurofeedback game could reward the user when their F3–P3 theta phase-locking value is high. This would require advanced signal processing and likely a multi-channel EEG setup, but it's a promising direction for a phase-based trainer. Even simpler, some systems train left-right **alpha asymmetry** (for depression/attention bias); for example, one PLOS Biology study in 2025 demonstrated that people could learn via feedback to **shift the asymmetry of their somatosensory alpha rhythms** (up-regulating alpha on one side vs the other) by focusing attention, mimicking the natural effect attention has on those oscillations <sup>36</sup>. They achieved bidirectional control: both focused attention and neurofeedback training were able to **increase or decrease alpha power asymmetry** on cue <sup>36</sup>. Although that study was examining pain perception, it provides a proof-of-concept that **people can volitionally modulate oscillation patterns (like lateralized alpha) relevant to attentional processing** when provided with real-time feedback.

**Neurofeedback training parameters:** Sessions are typically on the order of **20–30 minutes**, and training often extends over many sessions (e.g. 2–3 sessions per week for 5–10 weeks, totaling 20–40 sessions in clinical contexts). Changes in attention behavior usually emerge gradually, but can be long-lasting as they reflect learned self-regulation. A product aimed at healthy users might offer a shorter training regimen with clear goals (e.g. improve one's concentration span by X%, or reduce mind-wandering episodes). It could incorporate engaging game-like feedback (to keep motivation high during repetitive training). For example, a **focus training game** might display a character running when the user's beta:theta ratio is favorable, stopping when it isn't – turning abstract EEG metrics into a tangible challenge. Over time, the thresholds can be made more stringent as the user improves. Many modern neurofeedback systems also provide **progress tracking** (showing the user their trend in EEG metrics and attention test scores over sessions), which could be a feature of an attention phase trainer to quantify improvement.

One advantage of neurofeedback is that it can be combined with other modalities for synergy. For instance, **closed-loop tACS** could be paired with neurofeedback: the device could monitor EEG and only apply tACS when the user's attention dips, or conversely use neurofeedback to reinforce the brain after a tACS session (helping the user maintain the state that tACS induced). While such hybrid approaches are still experimental, they align with a product development mindset of maximizing efficacy by **integrating multiple techniques**.

In terms of evidence, neurofeedback's efficacy for attention is supported by numerous studies but also scrutinized. A *double-blind* trial by Eisenberg et al. (2021, for example) – not cited here – might show mixed results, whereas **open-label studies and meta-analyses** often find **significant improvements in inattention scores and cognitive performance** <sup>34</sup>. Importantly, **neurofeedback directly targets the user's own brain activity**, so any improvements are the user's achievement, potentially generalizing to everyday life without needing continued device use (after sufficient training). This contrasts with stimulation, where the effect may only last as long as you keep using the device. For a consumer attention trainer, offering neurofeedback means the product can aim for more sustainable, learned improvements in attention regulation.

## Design Considerations for an “Attention Phase Trainer” Tool

Drawing together the above findings, we can outline features and specifications for a non-invasive attention training device grounded in current neuroscience:

- **Multi-modal Stimulation Capability:** The device should support **tACS output** at least in the range of 1–50 Hz (covering delta up through gamma) with programmable waveforms (sine wave is standard, but options for turning stimulation on/off in bursts as used in some studies <sup>19</sup> would be useful). It should allow at least a **dual-electrode montage** (e.g. one frontal, one parietal) and ideally a **multi-channel montage** for flexibility (e.g. 4 electrodes so one could do two pairs or a ring). Critically, it must allow **phase offset control** between channels: e.g. the ability to set two channels 0° or 180° out-of-phase, or anywhere in between, with fine resolution. This is because research like Polanía's demands precise phase manipulation to recreate those effects <sup>5</sup>. Intensity control up to ~2 mA with safety monitoring is standard (most studies did 1–2 mA for 20 minutes safely <sup>6</sup>). The device would also include **EEG sensing** (dry or gel electrodes) so that it can read brain activity – necessary for neurofeedback and potentially for **adaptive stimulation** (e.g. tailoring frequency to the individual's peak as noted). Modern designs might integrate tACS and EEG on separate alternating current paths or in alternating time segments to avoid artifact, as simultaneous EEG

during tACS is challenging but methods exist to subtract the tACS artifact and recover EEG signals

37 .

- **Rhythmic Sensory Stimuli:** Include integrated **LED visual stimulators** (e.g. flickering light glasses or a screen interface) and **headphone or bone-conduction auditory output** for delivering rhythmic cues. These should be software-controlled with adjustable frequency and phase. For example, the user might follow a guided protocol where first a **10 Hz visual flicker** is presented for a few minutes (to entrain alpha), then switch to a **5 Hz auditory tone** (to train theta), etc., depending on the training program. The system could also offer **combined audio-visual entrainment** for stronger effects, but as research suggests, if both are used, one modality tends to dominate 32 . A strategy could be to alternate sessions (one session use visual, next session use tACS, etc.) or to use **sensory entrainment as an initial phase** to get the brain in rhythm, then turn it off and have tACS continue, or vice versa. The timing and sequence of multi-modal stimulation should be configurable, potentially using templates informed by research (e.g. a “gamma training” template might start with a 40 Hz flickering checkerboard for 1 minute to excite the visual cortex, then continue with 40 Hz tACS to frontal areas for 19 more minutes). Ensuring the **phase alignment** between audio and visual (if used together) is adjustable will allow exploration of optimum entrainment (for instance, sound clicks could be in phase with light flashes or offset by 180° to see which yields better attentional engagement).
- **Neurofeedback Interface:** The tool should provide a **real-time feedback display or game** driven by EEG metrics. Key EEG features to include are: **Theta/beta ratio (frontal midline)**, **alpha power (parietal or occipital)**, **beta or SMR power (central cortex)**, and possibly **frontal theta-parietal theta phase synchrony** if advanced signal processing is feasible. The software can offer various training “games” – e.g., one for TBR reduction (targeting improved sustained attention and reduced mind-wandering), one for alpha modulation (targeting either up-training for calm focus or down-training to reduce over-inhibition). For phase synchrony training, an interface might show two oscillating bars that become aligned when the user’s brain regions synchronize, and the user’s task is to “keep the bars aligned” via mental strategies. Incorporating the findings that **attention itself can modulate these rhythms** 36 , the system might include guided instructions like mindfulness or focusing exercises to help users find strategies to influence their EEG (e.g. “try to focus on a single point and clear your mind” may increase frontal midline theta, which some studies link to focused attention). Biofeedback could extend beyond EEG to include things like heart rate variability for a holistic training, but core to “phase trainer” is the EEG.
- **Protocol and Content:** Based on the literature, effective training likely requires **repeated practice over several weeks**. The product could suggest protocols such as: *Phase Training Level 1*: 3 sessions per week of 20 minutes tACS+entrainment + 10 minutes neurofeedback, for 4 weeks. Each session might start with a baseline attention test (like a 5-minute continuous performance task or go/no-go) to track progress. Then stimulation is applied while the user engages in an attention task (making it an **“attentional exercise” combined with entrainment** – this is important, as many studies delivered stimulation during task performance to effectively pair brain rhythm modulation with cognitive demands 38 19 ). For example, during 6 Hz frontoparietal tACS, the user could play a memory game or an attention-span game. After stimulation, the session could end with a neurofeedback cooldown, where the user, now without stimulation, tries to maintain the target brain state and sees their performance.



The inclusion of **specific reference protocols** is valuable: e.g., the device manual might cite “*Polanía 2012 protocol*” as one option (6 Hz, F3-P3, 0° phase, 15 min during working memory task) <sup>5</sup>, or “*Lee 2024 protocol*” (40 Hz, F4-P4, 2 mA, 20 min, repeated sessions) which was used to improve attention consistency in clinical subjects <sup>6</sup>. Having these as presets allows users or clinicians to choose evidence-based settings. Likewise, “*Alpha-FEF protocol*” (10 Hz intermittent, FEF targeting montage, 20 min) could be included, aiming to mimic top-down focus as per Misselhorn et al. 2024 <sup>19</sup>. For sensory entrainment, one preset could be “*Gamma flicker training*” (40 Hz visual flicker for 5 min + auditory 40 Hz for 5 min), drawn from findings that gamma RVS affects attention asymmetry <sup>28</sup>. Another could be “*Theta burst audio*” (a series of 5 Hz tone bursts, 2 seconds on, 2 seconds off, for 10 min) to practice rhythmic orienting of attention in time – inspired by research on rhythmic cueing. The device’s flexibility to implement these parameters with precision is crucial.

- **Monitoring and Adaptation:** Since individual responses vary, the tool should monitor EEG and performance data to **adapt the difficulty**. For neurofeedback, this means dynamically adjusting thresholds as the user improves (to keep the challenge optimal). For stimulation, adaptation could mean if a user isn’t responding to 10 Hz (no change in their EEG or behavior), the software might suggest trying 8 Hz or 12 Hz around that range, effectively finding a resonance (the importance of individual frequency tuning was highlighted when tACS had to target an individual’s preferred frequency for effect <sup>17</sup>). In a long-term training, one might even shift frequencies as the brain changes (maybe as attention improves, their dominant EEG frequency might increase – e.g., peak alpha might speed up slightly – and the protocol could follow that).
- **Outcome Assessment:** To validate improvement, the product should include **standardized attention assessments** before, during, and after the training regimen. These could be the Test of Variables of Attention (TOVA) as used by Lee et al. (which showed improvement in variability and impulsivity after gamma tACS <sup>15</sup>) or simpler continuous performance tests, go/no-go tasks (where reduction in omission and commission errors can be measured <sup>39</sup>). Even subjective questionnaires (mind-wandering frequency, ADHD symptom checklists) could be logged. By tying changes in these outcomes to the protocol parameters, the system could further refine what works best (for instance, noticing that users who did alpha entrainment show better distraction resistance, whereas those who did theta tACS show better working memory span, etc.).

In conclusion, the convergence of evidence from human studies indicates that **non-invasive rhythmic stimulation can be engineered to target the neural underpinnings of attention** in a precise way. Key frequencies – theta (~6 Hz) for executive control, alpha (~10 Hz) for sensory gating, beta/gamma (~15–40 Hz) for sustained alertness and processing – emerge as “dials” we can tune. Phase relationships (0° synchrony between frontal-parietal regions, or specific phase alignment between stimulus and brain oscillation) act as another dial for enhancing or suppressing connectivity <sup>3</sup>. A sophisticated attention trainer product would treat these like *training parameters*, much as a physical exercise machine lets one adjust resistance or speed. By leveraging tACS, visual/auditory entrainment, and neurofeedback in combination – all within safe, non-invasive limits – such a tool could provide **technical rigor in implementation (drawing on published protocols) and personalized adaptation**, thereby moving beyond generic brain-training into the realm of *evidence-based neuromodulation for cognitive enhancement*. With appropriate references to protocols <sup>3</sup> <sup>15</sup> <sup>28</sup> <sup>39</sup>, developers can trace the exact parameters that yielded benefits in studies (frequencies, electrode placements, session lengths) and incorporate those into the design. This ensures that the product is not just a gadget, but a translation of cutting-edge neuroscience into a practical attentional improvement regimen.

## Sources:

- Polanía et al. (2012), demonstrated causal impact of theta phase alignment on working memory <sup>3</sup> .
- Lee et al. (2024), 40 Hz tACS over right frontoparietal cortex improved attention consistency in clinical population <sup>15</sup> .
- Li et al. (2023), frequency-specific visual flicker (10 vs 40 Hz) altered visuospatial attention biases <sup>28</sup> .
- Cabral-Calderin & Henry (2025), rhythmic auditory vs. tACS entrainment – sensory rhythms dominated behavior, individual frequency tuning is key <sup>21</sup> <sup>32</sup> .
- Pahor & Jaušovec (2018), mixed effects of theta/gamma tACS on working memory, highlighting need for task-specific designs <sup>16</sup> .
- Hohn et al. (2025), neurofeedback enabled voluntary modulation of alpha oscillation asymmetry with attention, confirming trainability of oscillatory control <sup>36</sup> .
- Misselhorn et al. (2024), alpha-tACS over FEF mimicked top-down attention effects on visual processing <sup>11</sup> <sup>12</sup> .
- (Additional references embedded above provide further protocol details and results.)

---

<sup>1</sup> <sup>17</sup> <sup>18</sup> <sup>21</sup> <sup>22</sup> <sup>23</sup> <sup>31</sup> <sup>32</sup> <sup>37</sup> Sensory stimuli dominate over rhythmic electrical stimulation in modulating behavior | PLOS Biology

<https://journals.plos.org/plosbiology/article?id=10.1371/journal.pbio.3003180>

<sup>2</sup> <sup>16</sup> Frontiers | The Effects of Theta and Gamma tACS on Working Memory and Electrophysiology

<https://www.frontiersin.org/journals/human-neuroscience/articles/10.3389/fnhum.2017.00651/full>

<sup>3</sup> <sup>5</sup> <sup>7</sup> <sup>38</sup> The importance of timing in segregated theta phase-coupling for cognitive performance - PubMed

<https://pubmed.ncbi.nlm.nih.gov/22683259/>

<sup>4</sup> <sup>39</sup> Noninvasive Brain Stimulation for Improving Cognitive Deficits and Clinical Symptoms in Attention-Deficit/Hyperactivity Disorder: A Systematic Review and Meta-Analysis

<https://www.mdpi.com/2076-3425/14/12/1237>

<sup>6</sup> <sup>15</sup> Attention improvement to transcranial alternating current stimulation at gamma frequency over the right frontoparietal network: a preliminary report | Acta Neuropsychiatrica | Cambridge Core

<https://www.cambridge.org/core/journals/acta-neuropsychiatrica/article/abs/attention-improvement-to-transcranial-alternating-current-stimulation-at-gamma-frequency-over-the-right-frontoparietal-network-a-preliminary-report/933D462738B666CF48DB71DDD0290E5A>

<sup>8</sup> <sup>9</sup> <sup>10</sup> <sup>11</sup> <sup>12</sup> <sup>19</sup> <sup>20</sup> Transcranial Alternating Current Stimulation over Frontal Eye Fields Mimics Attentional Modulation of Visual Processing - PubMed

<https://pubmed.ncbi.nlm.nih.gov/38729759/>

<sup>13</sup> <sup>24</sup> <sup>25</sup> <sup>26</sup> <sup>27</sup> <sup>28</sup> <sup>29</sup> Frontiers | Modulation of rhythmic visual stimulation on left-right attentional asymmetry

<https://www.frontiersin.org/journals/neuroscience/articles/10.3389/fnins.2023.1156890/full>

<sup>14</sup> The Effect of Alpha tACS on the Temporal Resolution of Visual ...

<https://pmc.ncbi.nlm.nih.gov/articles/PMC7412991/>

<sup>30</sup> Alpha Power Increase After Transcranial Alternating Current ... - NIH

<https://pmc.ncbi.nlm.nih.gov/articles/PMC4464304/>

33 25-Month Follow-up of Double-Blind Randomized Controlled Trial

<https://pubmed.ncbi.nlm.nih.gov/36521694/>

34 a systematic review and meta-analysis - PMC - PubMed Central - NIH

<https://pmc.ncbi.nlm.nih.gov/articles/PMC9762555/>

35 Neurofeedback for ADHD: Time to Call It Quits? - Psychiatry Online

<https://psychiatryonline.org/doi/full/10.1176/appi.ajp.20220861>

36 Neurofeedback and attention modulate somatosensory alpha oscillations but not pain perception | PLOS Biology

<https://journals.plos.org/plosbiology/article?id=10.1371/journal.pbio.3002972>